

T: +44 (0) 203 930 8101 E: info@pcrbio.com Orders: sales@pcrbio.com

www.pcrbio.com

MATERIAL SAFETY DATA SHEET According to EC Directive 91/155 EC Date: 01/08/10

1. Identification of the substance/preparation and of the company/undertaking

1.1 Product Details:

Trade Names: PCRBIO Rapid Extract PCR Kit

Catalogue Numbers: PB10.24-08

PB10.24-40

Kit Components: 5x Rapid Extract Buffer A

10x Rapid Extract Buffer B 2x PCRBIO HS Tag Mix Red

1.2 Company Details

Manufacturer/Supplier:

PCR Biosystems Ltd. 1st Floor Aztec House 397-405 Archway Road London, N6 4ER United Kingdom

Telephone: +44 (0) 203 930 8101

Email: info@pcrbio.com

Further information obtainable from:

Technical Services Details as above

1.3 Emergency Details

24-hour contact in case of emergency:



T: +44 (0) 203 930 8101 E: info@pcrbio.com Orders: sales@pcrbio.com

www.pcrbio.com

PCR Biosystems Ltd. 1st Floor Aztec House 397-405 Archway Road London, N6 4ER United Kingdom

Telephone: +44 (0) 203 930 8101

2. Hazards identification

Product name Buffer B

Directive 1999/45/EC **IRRITANT**



R36/37/38 R42

GHS Directive 1272/2008/EC

DANGER



H315 - Skin Irrit. 2 H319 - Eye Irrit. 2 H334 - Resp. Sens. 1 H335 - STOT SE 3

2x PCRBIO HS

Tag Mix Red

Does not need labelling as

hazardous

Does not need labelling as

hazardous

According to Regulation (EC) No1272/2008 Skin Irritant (Category 2)

Eye Irritant (Category 2)

Respiratory Sensitivity (Category 1)

Specific Target Organ Toxicity (Single exposure) [Respiratory tract irritation] (Category 3)

According to European Directive 67/548/EEC as amended. Xi – Irritant

Label elements

Pictogram





Signal word

Warning Danger

Hazard statement(s)

H315 H319

H334

H335

Causes skin irritation.

Causes serious eye irritation.

May cause allergy or asthma symptoms or

breathing difficulties if inhaled.

May cause respiratory irritation.



T: +44 (0) 203 930 8101 E: info@pcrbio.com Orders: sales@pcrbio.com

Precautionary statement(s)
P261S Avoid breathing dust

P280sh Wear protective gloves/eye protection
P302+352 IF ON SKIN: Wash with plenty of soap and water
P305+351+338 IF IN EYES: Rinse cautiously with water for

several minutes. Remove contact lenses, if present and easy to

do. Continue rinsing.

Hazard symbol(s)

Xi Sensitising

R-phrase(s)

R36/37/38 Irritating to eyes, respiratory system and skin.

R42 May cause sensitisation by inhalation.

S-phrase(s)

S36/37 Wear suitable protective clothing and gloves

Other hazards None

3. Composition/information on ingredients

· Chemical characterization: Mixtures

· Description: Mixture of the substances listed below with nonhazardous additions.

· Dangerous components:

Glycerol CAS: 56-81-5 EINECS: 200-289-5 Proteinase K

EC-No 254-457-8, 5-10%, Classification H315;H319;H334;H335

· Additional information For the wording of the listed risk phrases refer to section 16.

4. First aid measures

- · Description of first aid measures
- · General information No special measures required.
- · After inhalation Supply fresh air; consult doctor in case of complaints.
- · After skin contact Generally the product does not irritate the skin.
- · After eye contact Rinse opened eye for several minutes under running water.
- After swallowing Induce vomiting and call for medical help.
- · Information for doctor
- · Most important symptoms and effects, both acute and delayed No further relevant information available
- · Indication of any immediate medical attention and special treatment needed No further relevant information available.

5. Firefighting measures

- Extinguishing media
- · Suitable extinguishing agents

CO2, extinguishing powder or water spray. Fight larger fires with water spray or alcohol resistant foam.

- · Special hazards arising from the substance or mixture No further relevant information available.
- Advice for firefighters
- · Protective equipment: No special measures required.



T: +44 (0) 203 930 8101 E: info@pcrbio.com Orders: sales@pcrbio.com

www.pcrbio.com

6. Accidental release measures

- · Personal precautions, protective equipment and emergency procedures Not required.
- · Environmental precautions:

Dilute with plenty of water.

Do not allow to enter sewers/ surface or ground water.

Methods and material for containment and cleaning up:

Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).

Dispose contaminated material as waste according to item 13.

· Reference to other sections

No dangerous substances are released.

See Section 7 for information on safe handling

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

7. Handling and storage

- · Handling
- · Precautions for safe handling No special measures required.
- · Information about protection against explosions and fires: No special measures required.
- · Conditions for safe storage, including any incompatibilities
- Storage
- · Requirements to be met by storerooms and receptacles: According to product specification
- · Information about storage in one common storage facility: Not required.
- Further information about storage conditions: None.
- · Specific end use(s) No further relevant information available.

8. Exposure controls/personal protection

- · Additional information about design of technical systems: No further data; see item 7.
- · Control parameters
- · Components with limit values that require monitoring at the workplace:

56-81-5 glycerol

WEL () Long-term value: 10 mg/m³

- · Additional information: The lists that were valid during the creation were used as basis.
- · Exposure controls
- · Personal protective equipment
- · General protective and hygienic measures

The usual precautionary measures should be adhered to when handling chemicals.

- · Breathing equipment: Not required.
- · Protection of hands: Protective gloves.
- · Material of gloves Synthetic gloves
- Penetration time of glove material

The exact break trough time has to be found out by the manufacturer of the protective gloves and has to be observed.

· Eye protection: Safety glasses

9. Physical and chemical properties

- · Information on basic physical and chemical properties
- · General Information



simplifying research

PCR Biosystems Ltd Aztec House 397-405 Archway Road London, N6 4ER

T: +44 (0) 203 930 8101 E: info@pcrbio.com Orders: sales@pcrbio.com

www.pcrbio.com

Appearance:

Form: Fluid
Colour: Colourless
· Odour: Odourless

· Odour threshold: Not determined.

· **pH-value:** Not determined.

· Change in condition

Melting point/Melting range: undetermined

Boiling point/Boiling range: 100°C

Flash point: > 55°C

· Flammability (solid, gaseous) Not applicable.

· Ignition temperature: 270°C

· Decomposition temperature: Not determined.

· **Self igniting:** Product is not selfigniting.

· Danger of explosion: Product does not present an explosion hazard.

• Explosion limits: Lower: 1.8 Vol % Upper: 63.0 Vol %

· Vapour pressure at 20°C: 23 hPa

· Density: Not determined

Relative density Not determined.

· Vapour density Not determined.

• Evaporation rate Not determined.

· Solubility in / Miscibility with

Water: Fully miscible

· **Segregation coefficient (n-octanol/water):** Not determined.

· Viscosity:

dynamic: Not determined. **kinematic:** Not determined.

· Solvent content:

Organic solvents: 12.0 %

Water: 71.3 %

Solids content: 16.7 %

· Other information No further relevant information available.

10. Stability and reactivity

Reactivity

· Chemical stability

- · Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
- · Possibility of hazardous reactions No dangerous reactions known
- · Conditions to avoid No further relevant information available.
- · Incompatible materials: No further relevant information available.
- · Hazardous decomposition products: No dangerous decomposition products known

11. Toxicological information

- · Information on toxicological effects
- · Acute toxicity:
- · LD/LC50 values that are relevant for classification:

67-68-5 dimethyl sulfoxide

Oral LD50 14500 mg/kg (rat)



T: +44 (0) 203 930 8101 E: info@pcrbio.com Orders: sales@pcrbio.com

www.pcrbio.com

· Primary irritant effect:

· on the skin: No irritant effect.

- · on the eye: No irritating effect.
- · Sensitization: No sensitizing effects known.
- · Additional toxicological information:

The product is not subject to classification according to the calculation method of the General EU Classification Guidelines for Preparations as issued in the latest version:

When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.

12. Ecological information

- · Toxicity
- · Acquatic toxicity: No further relevant information available.
- · Persistence and degradability No further relevant information available.
- · Behaviour in environmental systems:
- · Bioaccumulative potential No further relevant information available.
- · Mobility in soil No further relevant information available.
- Additional ecological information:
- · General notes:

Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water.

Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

- · Results of PBT and vPvB assessment
- · **PBT:** Not applicable.
- · **vPvB**: Not applicable.
- · Other adverse effects No further relevant information available.

13. Disposal considerations

- · Waste treatment methods
- · Recommendation

Must not be disposed of together with household garbage. Do not allow product to reach sewage system.

- · Uncleaned packagings:
- · Recommendation: Disposal must be made according to official regulations.
- · Recommended cleansing agent: Water, if necessary with cleansing agents.

14. Transport information

- Maritime transport IMDG:
- · Marine pollutant: No
- · Special precautions for user Not applicable.
- · Transport/Additional information: Not dangerous according to the above specifications.
- · Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code Not applicable.

15. Regulatory information

- · Safety, health and environmental regulations/legislation specific for the substance or mixture
- · National regulations
- · Technical instructions (air):

Class Share in %

NK 10-20

- · Water hazard class: Water hazard class 1 (Self-assessment): slightly hazardous for water.
- · Chemical safety assessment: A Chemical Safety Assessment has not been carried out.



T: +44 (0) 203 930 8101 E: info@pcrbio.com Orders: sales@pcrbio.com

www.pcrbio.com

16. Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.